


# A large and versatile range of tips interacting in harmony with the handpiece and the device to deliver optimum performance

Newtron® tips are conceived to meet all clinical needs, thanks to exclusive designs, alloys and coatings that respect the surfaces treated: enamel, crown, implant.





# scaling



F00934

Supra- and sub-gingival scaling  
N° 1, N° 1S, N° 10X, H3 tips,  
4 autoclavable dynamometric  
wrenches



1

F00246

2

F00247



3

F00248

10P

F00253



10Z

F00254

## Supra-gingival scaling



### Universal tip

Simple cases: gross supra-gingival scaling.

Tangential orientation to the surface.  
To-and-fro sweeping to "detach" the tartar whilst respecting the enamel.



### Voluminous calculus

Removal of significant supra-gingival deposits.

Apply the flat part to the tooth surfaces.



### Stains

Removal of marks and stains (tobacco, tea, coffee, etc.).

Apply the rounded extremity of the tip to the surface to be treated.

## Sub-gingival scaling and probing



### Shallow pockets

Scaling of pockets less than 2-3mm deep.




### Medium pockets

Scaling of medium pockets (< 4mm).  
Removal of biofilm and soft deposits, while evaluating the depth of the pockets using the marks every 3mm.  
Efficient for maintenance treatment in patients with good dental hygiene.



# hygiene



F00935

N° 1, N° 1S, N° 10Z, TK1-1S tips,  
4 autoclavable dynamometric  
wrenches



## Supra- and sub-gingival scaling



### Slim tip

Interproximal spaces scaling. Finer and longer than tip No.1, it is also powerful and robust.

## Supra-gingival scaling and interproximal spaces



### Interproximal spaces

Its anatomical shape allow fast and efficient procedure.



1S

F00245



10X

F00359

# periofine

## Smooth biofilm elimination



### Dental plaque and sub-gingival small deposits removal

Oriented tangentially: its shape adapts to the anatomy of the tooth for a painless and easy access.



### Interproximal scaling of narrow areas

Left-oriented for an easy access to premolars and molars.



### Interproximal scaling of narrow areas

Right-oriented for debridement and cleaning of medium pockets.



PFU

F02170



PFL


F02171




PFR

F02172


# periodontics



F00936 N° 1S, H3, H4L, H4R tips,  
4 autoclavable dynamometric  
wrenches



## Periodontal debridement



H3

### Initial periodontics, anterior sector

Treatment of the incisor-canine block.

The guide edge is oriented parallel to the pocket. The H3 tip is descended into the periodontal pocket without risk of injury to the ligament. The cavitation will lift the debris out.



H4L

### Periodontics for the premolar and molar sectors, left-oriented

First instrument in the sequence for treating all the surfaces and the furcations.

- Maxillary: buccal and distal surfaces of sector 2, pivots at 13, then the buccal and mesial surfaces of sector 1.
- Mandibular: buccal and distal surfaces of sector 4, pivots at 43, then lingual and mesial surfaces of sector 3.



H4R

### Periodontics for the premolar and molar sectors, right-oriented

Second instrument in the sequence.

- Maxillary: palatine and mesial surfaces of sector 2, pivots at 13, then buccal and distal surfaces of sector 1
- Mandibular: lingual and mesial surfaces of sector 4, pivots at 43, then buccal and distal surfaces of sector 3.

## Root planing



H1

### Anterior tooth root planing, diamond-coated tip 30 µm

- Diamond-coated mini-tip for simple cases in the cervical area.
- Also effective for the withdrawal of granulation tissue.


This tip should be used without pressure and above the epithelial attachment because it is abrasive.



H2L

### Root planing of the premolar and molar sectors, left-oriented, diamond-coated tip 30µm

Diamond-coated micro-probe for the treatment of furcations and narrow spaces.



H2R

### Root planing of the premolar and molar sectors, right-oriented, diamond-coated tip 30µm

Diamond-coated micro-probe for the treatment of furcations and narrow spaces.



H1  
F00366



H2L  
F00367



H2R  
F00368



H3  
F00369




H4L  
F00114



H4R  
F00115



# perio maintenance BDR



F00737

TK1-1S, TK1-1L, TK2-1L, TK2-1R  
tips, 4 autoclavable dynamometric  
wrenches



# perio precision



F00939

P2L, P2R, TK1-1S tips,  
3 autoclavable dynamometric  
wrenches



## Biofilm disruption



### Short probe

Graduated every 3mm, for examining shallow and medium pockets (< 4mm) and for the maintenance of simple cases.



### Long probe

Examination and maintenance of medium to deep pockets (> 4mm).  
Diagnosis aid during the debridement and irrigation of pockets.

The TK1 probe tips are used without pressure following the contour of the pockets and skimming over the root surface.



### Maintenance of the premolar and molar sectors, left-oriented


Maintenance of moderate to deep pockets and furcations.  
Equivalent to the Nabers probe.



### Maintenance of the premolar and molar sectors, right-oriented

Complementary to the TK2-1L tip for the maintenance of moderate to deep pockets and furcations.  
Equivalent to the Nabers probe.

## Periodontal maintenance



### Debridement of the premolar and molar sectors, left-oriented

Round micro-tip recommended for periodontal debridement in the presence of a fine peridontium and in narrow areas.

- Maxillary: buccal and distal surfaces of sector 2, pivots at 13, then the palatine and mesial surfaces of sector 1.
- Mandibular: buccal and distal surfaces of sector 4, pivots at 43, then lingual and mesial surfaces of sector 3.




### Debridement of the premolar and molar sectors, right-oriented

Second instrument in the sequence, after the P2L tip.

The double bend makes it possible to treat areas that are difficult to access (inter-radicular spaces, deep pockets).

- Maxillary: buccal and mesial surfaces of sector 2, pivots at 13, then buccal and distal surfaces of sector 1.
- Mandibular: lingual and mesial surfaces of sector 4, pivots at 43, then buccal and distal surfaces of sector 3.

The P2 tips can also be used to remove small amounts of excess cement when bonding fixed prosthesis.



→ TK2-1L / P2L  
→ TK2-1R / P2R



TK1-1S  
F01001



TK1-1L  
F01004



TK2-1L  
F02162



TK2-1R  
F02161



P2L  
F00090



P2R  
F00091

# periosoft



## Implant and prosthesis prevention



### Hygiene of anterior sector

Plastic micro-tip with universal curette shape for the treatment of the incisor/canine groups.

- Removal of the biofilm and low adherence deposits without scratching the prosthetic surfaces.
- Polishing the sulcus or grooves of natural teeth.



### Hygiene of premolar and molar sectors, left-oriented

Plastic micro-tip with 13-14 curette shape for the removal of biofilm and low adherence deposits for the treatment of the posterior groups.

- Maintenance for the screws and abutment of the implant.
- Scaling of prosthesis.



### Hygiene of premolar and molar sectors, right-oriented

Plastic micro-tip with 13-14 curette shape for the removal of biofilm and low adherence deposits for the treatment of the posterior groups.

*The new material for these tips makes it possible to clean and debride faster, and gives better breakage resistance.  
Max. Power = 3 (start of green mode).*


# implantprotect pure titanium

*Pure titanium tips to preserve implant surfaces.*

F02120



*IP1, IP2L, IP2R, IP3L, IP3R tips,  
autoclavable metal support and  
universal wrench*



## Treatment of peri-implantitis and maintenance




### Debridement of the implant abutment and wide threads

Pure titanium tip with a wider extremity for implant abutment cleaning and large thread debridement.



### Debridement of medium implant threads, left-oriented

Pure titanium tip with a similar shape to P2L tip for the debridement of medium implant threads. The bend of the tip allows movement around the entire implant for total decontamination.



### Debridement of medium implant threads, right-oriented

Pure titanium tip with a similar shape to P2R for the debridement of medium-sized implant threads. The approach may be non-surgical or open flap.



### Debridement of narrow implant threads, left-oriented

Pure titanium tip with a pointed extremity suitable to reach narrow implant threads. All types of implants can be treated with these different tip sizes.



### Debridement of narrow implant threads, right-oriented

Pure titanium tip with a pointed extremity suitable to reach the inner-most parts of narrow implant threads.

*The black ring on these tips indicates their exclusive use on titanium.  
Max. Power = 5 (green)*



# endosuccess canal access prep

The micro-blades are less aggressive than diamond and their coating makes these tips very durable.

F88180



CAP1, CAP2, CAP3 tips, autoclavable metal support and universal wrench



## Canal access preparation



CAP1

### Micro-blade tip length 12mm, taper 6%

Active lateral part for:

- Finishing walls and polishing.
- Removing temporary cement and dentinal residues.
- Removing dentin overhangs.

Non-active end to prevent the risk of perforating the pulp chamber floor.



CAP2

### Micro-blade tip, length 9mm, taper 5%

Active lateral part and extremity used by sweeping method to remove dentine bridges.

- Location of the MB2 (2nd mesiobuccal canal) and search for hidden canals.
- Preparation of the pulp chamber.
- Removal of the dentine layer which may hide the access to the MB2 canal.



CAP3

### Micro-blade tip, length 8mm, taper 6%

The CAP3 tip has a very pointed extremity indicated for:

- Locating and opening the calcified canals.
- Fragmenting calcifications or pulp stones in the pulp chamber.
- Loosening fiber posts.
- Locating accessory canals.

Due to its very sharp point, the CAP3 tip must be handled with care (visual aids recommended).



ET18D

### Diamond-coated steel tip 76µm, length 18mm, taper 5%

- Finishing the access cavity.
- Removing dentine overhangs, calcifications and filling materials.



ETBD

### Diamond-coated ball tip, length 20mm, taper 5%

Searching for canals and locating calcified canals.

CAP1  
F88181CAP2  
F88182CAP3  
F88183ET18D  
F88017ETBD  
F88020

# irrisafe



IRRI20,25

## Irrigation



IRRISAFE

### Passive ultrasonic irrigation (PUI) files of different lengths and diameters

Irrisafe™ safely\* eliminates the smear layer, dentine debris and bacteria from the root canal. Its blunt tip prevents any risk of perforating the apex or the canal walls.

Irrigation once the root canal has been prepared.

- 20ml of irrigant (NaOCl) are injected into the canal.
- Irrisafe™ is inserted 2mm short of the working length and activated by performing withdrawal movements to flush the debris and the smear layer upwards.
- Repeated 3x 1 minute in each canal.



K FILES

### Files of different lengths and diameters, taper 2%

Irrigation, withdrawal of calcified dentine and gutta percha, and withdrawal of broken instruments.

For irrigation ultrasonic files are used with a disinfectant solution. To provide a final decontamination, use sodium hypochlorite until the smear layer is removed.

K files are very sharp instruments and should be handled with precision. However they are flexible and can therefore be pre-bent.

K10, 15, 25, 30  
FILES

\* Van Der Sluis L.W.M. Passive ultrasonic irrigation of the root canal: a review of the literature. Int. Endodont. J. 2007; 40; 4: 415-428



# endosuccess retreatment

F00737



ET18D, ET20, ET25, ET25S, ETBD,  
ETPR tips, autoclavable metal support  
and universal wrench



# endodontics

F00732  
endo  
one

Endodontic treatments  
CAP1, CAP2, CAP3, ET25, ETPR  
tips, 4 Irrisafe 25-21 mm blister,  
autoclavable metal support  
and universal wrench



## Canal Retreatment



ET20

### Retreatment tip, length 20mm, taper 6%

Used in the 1<sup>st</sup> coronal third:

- Extraction of filling material, silver points, broken instruments.
- Removal of debris and the smear layer.



ET20D

### Diamond-coated retreatment tip, 30 µm, length 20mm, taper 5%

Used in the 1<sup>st</sup> coronal third to remove very hard materials by brushing the walls. The diamond coating of the ET20D tip increases the cutting and lateral abrasion effect.



ET25

### Titanium-Niobium tip, length 20mm, taper 3%

Retreatment in the middle and apical thirds and the extraction of broken instruments. The Titanium-Niobium alloy of the ET25 range allows perfect transmission of the ultrasonic vibrations and tip flexibility\*.



ET25S

### Short Titanium-Niobium tip, length 15mm, taper 4%

Retreatment in the coronal third and the isthmuses.

## Retreatment and obturation



ET40

### Long retreatment tip, 40mm, taper 4%


Rapid removal of broken instruments in the middle third of wide, straight canals.



ET40D

### Long retreatment tip, 40mm, diamond-coated 30 µm, taper 4%

Retreatment of very hard material in the middle third.



ET25L

### Long Titanium-Niobium tip, 25mm, taper 3%

Retreatment in the apical third and long, straight canals.


ET25 tips can be pre-formed for the treatment of curved canals.



SO4

### Fine condenser, length 40mm, taper 4%

Lateral condensation of gutta percha by heating effect. It is used dry, without irrigation.

ET20  
F88011ET20D  
F88013ET25  
F88018ET25S  
F88021ET40  
F88012ET40D  
F88014ET25L  
F88022SO4  
F88009

\* E.W. Collings Applied superconductivity, metallurgy and physics of titanium alloys 1985

# endosuccess apical surgery

F00069



AS3D, AS6D, AS9D, ASLD, ASRD  
tips, autoclavable metal support and  
universal wrench



## Apical surgery



**Diamond-coated universal tip 30μm, length 3mm, taper 9%**

Apical surgery of anterior teeth.  
It should be used without pressure, at the lowest possible effective power.



**Diamond-coated tip 30μm, length 6mm, taper 9%**

Second instrument in the sequence, used to obtain a preparation length of 5mm at least.



**Diamond-coated tip 30μm, length 9mm, taper 8%**

Used for complex cases and for the preparation of the root canal up to the coronal third.  
The diamond coating is only present on the extremity of the instrument not to over-prepare the canal.

*The AS9D tip should first be introduced into the canal and oriented in the root axis before being activated to prevent the creation of a «false route».*



**Right-oriented tip, diamond-coated 30μm, length 3mm, taper 10%**

Apical surgery of premolars and molars.



**Left-oriented tip, diamond-coated 30μm, length 3mm, taper 10%**

Apical surgery of premolars and molars.

*It should be used with very light pressure.*



AS3D  
F00065



AS6D  
F00079



AS9D  
F00067




ASRD  
F00081

ASLD  
F00080

# endosurgery

## Retro surgery



**Retro surgery tip angled at 70°, diamond-coated 30μm, length 5mm, taper 9%**

Treatment of posterior areas, in canals that are difficult to access or roots with specific orientations.



**Universal retro surgery tip, diamond-coated 30μm, length 5mm, taper 7%**

Preparation of canals in anterior teeth.  
The micro-retro tips make minimum treatment possible providing fast healing.



**Left-oriented retro surgery tip, diamond-coated 30μm, length 5mm, taper 7%**

Preparation of premolar and molar canals.



**Right-oriented retro surgery tip, diamond-coated 30μm, length 5mm, taper 7%**

Preparation of premolar and molar canals.



S12-70D  
F00118



P14D  
F00106



P15LD  
F00107

P15RD  
F00108



# perfect margin rounded

F00738



PM1, PM2, PM3, PM4 tips,  
autoclavable metal support and  
universal wrench



## Prosthetic finishing with chamfered shape



PM1

### Preparation, rounded edge, diamond-coated tip 76 µm

First instrument of the ultrasonic sequence, following the rotary phase. Intraculcular dentin preparation and positioning of finishing line.



PM2

### Finishing, rounded edge, diamond-coated tip 46 µm

Correction of irregularities in the finish line and start of polishing. Its diamond coating, less dense than on the PM1, makes it possible to obtain a cutting edge finish.



PM3

### Polishing, rounded edge, smooth

This entirely smooth instrument is last in the finishing sequence, improving the condition of the surface at the cervical limit before impression taking.



PM4

### Corono-radicular preparation, conical, diamond-coated 46 µm

After the rotating phase, the PM4 tip is used to:

- Prepare the upper 1/3 of canal chamber.
- Shape anatomically the connection cone.
- Clean the root walls.
- Smooth the entry cones for the anatomical posts.

## Loosening and condensation

5AE

### Loosening of root canal posts with spray

Apply the 5AE tip on the lingual or palatine surface and the buccal surface, before finishing with the occlusal surface. Use the flat extremity of the instrument held firmly against the tooth.

C20

### Condensation, Piezocem

For inlays or onlays on posterior teeth.

ETPR

### Loosening tip (post removal)

The ETPR tip has profiled and concave shape. It provides greater efficacy on the posterior teeth.



5AE  
F00249



C20  
F00113



ETPR  
F88019