

Dental Equipment

pages 8-15 Cleo II Treatment Centres



Cleo II

Advanced technology, improved ergonomics and the Operator-friendly features of the Cleo II ensure superior control and precision for the dental team. The flexibility of the Cleo II range ensures that any configuration will equip a dental surgery with the ultimate, complete surgery system. The unique folding and extending chair provides unrivalled comfort and access for both Operator and Patient.

CLEO II 'Surgery System'
Air or Electric versions





CLEO II 'Flexible System' (A) or (E)



Mobile Cart
Air or Electric versions

Cabinet Module Air or **Electric** versions

Patient-friendly, synchronised chair:

Unique folding-extending chair occupies an unusually small footprint and provides unrivalled comfort, natural support and enhanced access for both Operator and Patient.



uniquely ideal:



Equipped for an array of surgery procedures and working methods in complete comfort.

Height-adjustable, extendable, rotating Assistant's console Easy-clean porcelain spittoon with cup-filler and bowl-flush Rotating, chair-mounted cuspidor (90°) 'E' version key features ~ Rotating cuspidor unit LED digital display increases access for an Assistant working closely Preset micro motor Torque control with Auto-Reverse over the Patient from the 3 o'clock position First Priority handpiece selection Operator's console is height-adjustable and rotates: supreme versatility and the pneumatic brake enables precise positioning for ideal working comfort. Removable rubber mat **Control and precision** are the prerequisites of the Cleo II. The Operator's console rotates smoothly on a compact, base-mount enabling an Operator easy access between the 9 and 1 o'clock positions. Hygiene-conscious at every point, the Cleo II is designed to perform in the most demanding dental surgery environments.

Versatility for working comfort:

From consultation through treatment, Cleo II provides supreme access and working comfort for the Operator and Assistant.



Folding legrest and detachable, foldaway armrest for perfect Patient entry and exit

Clean and preparation during consultation improves efficiency




Maximised workflow and a greater variety of operation are the key acheivements of the Cleo II treatment centre. An inclined Cleo II chair is a very gentle welcome to any Patient entering the surgery and enables the most pleasant consultation position possible.



complete, hygienic surgery system:

Cleo II



Cleo | 'Flexible' treatment centres :

Mobile configurations to suit any Practice and the flexibility of Electric and Air versions to fulfil individual preferences.

Tubes are easily detached

CLEO II 'FLEXIBLE SYSTEM'

Rotating cuspidor (90°)

Foldaway, detachable armrest



Ciso II chair key features ~
Folding and extending legrest
Twin-articulating headrest
Foldway, detachabele armest
Synchronised chair movements
Lumbar support
Seamless upholstery

Folding and extending legrest safety
feature ensures that the legrest retracts
if contact is made with any object

Integrated foot control stalks for chair movements including two programmable presets, Last-Position memory

and Auto-Return

Ultra-slim backrest improves Operator's working comfort

The telescopic base is a superior, elegant design solution for easy hygiene maintenance – smooth, white surfaces, with integrated foot controls. Supremely comfortable chair movements, smooth, near-silent and gentle start/stop.

technically advantaged:

Cleo II

Super ergonomics for Dentist and Patient.

'E' version control panel includes instrument, chair, spittoon and light controls



'A' version control panel includes chair and spittoon controls





Mobile Cart and Cabinet Module 'A' and 'E' versions



Increased range of synchronised chair movements



Assistant's console with chair, spittoon and light¹ controls

Control all chair movements by hand or foot: by hand from both the Operator's and Assistant's console; by foot using the stalks built-in to the chair base. Spittoon functions and ON/OFF light control $\ensuremath{^{1}}$ are also activated from both the Operator's and Assistant's console.



820S LIGHT - Unit-mounted or 700 SERIES LIGHT - Unit- or ceiling-mounted

Triple-axis head movement White natural beam 5,000 Kelvin Touchless sensor ON/OFF control Composite Curing Safe Mode

4,200 Kelvin Touchless sensor ON/OFF control Composite Curing Safe Mode Retractable Patient mirror

Operator's foot control for instruments